

INSPIRE* RECOMBINANT HUMANTM TYPE-I COLLAGEN WOUND DRESSING, A NOVEL SCAFFOLD FOR CUTANEOUS WOUND HEALING.

Shani Shilo¹, Sigal Roth¹, Or Dgany¹, Olga Mizrahi² B.Med.Sc., Dima Sheyn² MSc., Gadi Pelled² PhD, DMD, Dan Gazit² PhD, DMD, Hagit Amitai¹, Oded Shoseyov^{1,3}

¹CollPlant Ltd. 3 Sapir St. Weizmann Science Park P.O.B 4132 Ness-Ziona 74140, Israel

²Skeletal Biotech Laboratory, Hebrew University of Jerusalem – Hadassah Medical Campus, Jerusalem, 91120, Israel

³The Robert H. Smith Institute of Plant Science and Genetics and the Minerva Otto Warburg Center for Agricultural Biotechnology. The Robert H. Smith Faculty of Agricultural, Food and Environmental Quality Sciences, the Hebrew University of Jerusalem P.O.B. 12 Rehovot 76100, Israel.

Type-I Collagen scaffolds are thought to be the best "man made" scaffolds for the medical device market, providing an optimal mimic environment for cell proliferation and tissue repair. To date, the collagen used is of animal origin, entailing several risks including pathogens and allergic reactions. The purpose of this study was to produce a safe and efficient scaffold based on novel human recombinant type-I collagen. Using genetic engineering we have successfully expressed five human genes in tobacco plants to produce fully functional triple helical recombinant human type-I collagen, **Collage-rhTM**, identical to human collagen. A 100% human recombinant collagen type I wound dressing sheet, "**Inspire-rhTM**" wound dressing sheet was fabricated by fibrillogenesis and freeze drying. Thermal dehydration treatment was applied to physically cross-link the collagen monomers in order to increase chemical and mechanical stability. A highly porous cohesive structure was obtained, capable of absorbing large quantities of fluids, up to 40 times its own weight, in a few seconds. This characteristic is particularly important and highly desirable for wound dressing products, due to the need to absorb exudates. In tissue culture the bio-functionality of "**Inspire-rhTM**" wound dressing sheet was studied using key cell models that participate in the wound healing process, human dermal fibroblasts, human endothelial cells and human bone marrow-derived mesenchymal stem cells. Proliferation of cells on **Inspire-rhTM** scaffolds was compared to commercially available collagen containing products for cutaneous wound healing. **Inspire-rhTM** sheets showed superior biological performance in supporting human dermal fibroblasts, human endothelial cells and human mesenchymal stem cell proliferation, therefore expected to jump start tissue repair and accelerate wound healing. **Collage-rhTM** can be used to produce additional types of safe and efficient scaffolds for other applications such as orthopedic bone devices, leading to a new era in regenerative medicine.

***Inspire-rhTM** was later changed to **VergenixTM**